

1 mL, Model 1001 Saltline (SAL), Instrument Syringe

Specification Sheet (Part/REF # 208335)



The Bubble Free Prime syringes are compatible with all Microlab 600 applications. The SaltLine syringes were specially developed to withstand harsh working conditions in demanding applications over a long product lifetime.

Product Specifications

Needle Sold Separately	Yes
Volume	1 mL
Syringe Type	Gastight, Saltline
Syringe Series	1000
Needle	No Needle Available
Plunger Style	Saltline (SAL) X-Style
Fluid Path	Borosilicate glass, UHMWPE, PTFE, CTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack
Plunger Tip Material	UHMWPE
Thread at Termination	5/16-28

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/208335>

Spec. Version B