

# 10 mL, Model 1010 Saltline (SAL), Instrument Syringe

Specification Sheet (Part/REF # 208337)

The Bubble Free Prime syringes are compatible with all Microlab 600 applications. The SaltLine syringes were specially developed to withstand harsh working conditions in demanding applications over a long product lifetime.

## Product Specifications

<b>Needle Sold Separately</b>	Yes
<b>Volume</b>	10 mL
<b>Syringe Type</b>	Gastight, Saltline
<b>Syringe Series</b>	1000
<b>Needle</b>	No Needle Available
<b>Plunger Style</b>	Saltline (SAL) X-Style
<b>Fluid Path</b>	Borosilicate glass, UHMWPE, PTFE, CTFE
<b>Autoclavable</b>	No
<b>Gas Sterilizable</b>	Yes
<b>Pressure</b>	200 psig (13.8 bar)
<b>Maximum Temperature</b>	50°C / 122°F
<b>Minimum Temperature</b>	20°C / 68°F
<b>Sales Quantity</b>	1 syringe per pack
<b>Plunger Tip Material</b>	UHMWPE
<b>Thread at Termination</b>	5/16-28

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/208337>

**Spec. Version B**