

5 mL, Model 1005 Saltline (SAL), Instrument Syringe

Specification Sheet (Part/REF # 208336)



The Bubble Free Prime syringes are compatible with all Microlab 600 applications. The SaltLine syringes were specially developed to withstand harsh working conditions in demanding applications over a long product lifetime.

Product Specifications

| | |
|-------------------------------|--|
| Needle Sold Separately | Yes |
| Volume | 5 mL |
| Syringe Type | Gastight, Saltline |
| Syringe Series | 1000 |
| Needle | No Needle Available |
| Plunger Style | Saltline (SAL) X-Style |
| Fluid Path | Borosilicate glass, UHMWPE, PTFE, CTFE |
| Autoclavable | No |
| Gas Sterilizable | Yes |
| Pressure | 200 psig (13.8 bar) |
| Maximum Temperature | 50°C / 122°F |
| Minimum Temperature | 20°C / 68°F |
| Sales Quantity | 1 syringe per pack |
| Plunger Tip Material | UHMWPE |
| Thread at Termination | 5/16-28 |

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/208336>

Spec. Version B