

12.5 µL, Model 1702.5 TLLX SYR, Instrument Syringe

Specification Sheet (Part/REF # 5495-10)



Product Specifications

Volume	12.5 µL
Syringe Type	Gastight Half-Height
Syringe Series	1700
Termination	PTFE Luer Lock (TLL)
Plunger Style	X-Style
Fluid Path	Borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.029 inches (0.729 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/5495-10>

Spec. Version B