

2.5 mL, Model 1005.5 TLL SYR, Long Life Instrument Syringe

Specification Sheet (Part/REF # 8250-45)



Product Specifications

Volume	2.5 mL
Syringe Type	Longlife Gastight Half-Height
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Plunger Style	X-Style
Fluid Path	Borosilicate glass, UHMWPE, PFA
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.406 inches (10.3 mm)
Barrel Outer Diameter	0.530 inches (13.5 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/8250-45>

Spec. Version B