

10 µL, Model 1701 RN Agilent SYR, Small Removable NDL, 23s-26s ga, 1.71 in, point style AS

Specification Sheet (Part/REF # 80089)



Product Specifications

Volume	10 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	23s-26s gauge
Point Style	AS: Conical, Autosampler
Needle Length	1.71 inch (43.4 mm)
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.018 inches (0.460 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Needle Inner Diameter	0.0133 - 0.005 inches (0.64 - 0.13 mm)
Needle Outer Diameter	0.0283 - 0.0187 inches (0.72 - 0.47 mm)
Std Needle Dead Volume	0.46 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/80089>

Spec. Version B