

10 µL, Model 701 N Agilent SYR, Cemented NDL, 26s ga, 1.71 in, point style 2

Specification Sheet (Part/REF # 80399)



Hamilton syringes designed for Agilent GC Autosamplers 7693A, 7673, 7683, and 6850 ALS feature a specially designed plunger button.

Product Specifications

Volume	10 µL
Syringe Type	Microliter
Syringe Series	700
Termination	Cemented Needle (N)
Needle	Cemented Needle
Gauge	26s gauge
Point Style	2: Curved, Beveled
Needle Length	1.71 inch (43.4 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	No
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.019 inches (0.485 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Needle Inner Diameter	0.005 inches (0.13 mm)
Needle Outer Diameter	0.0187 inches (0.47 mm)
Std Needle Dead Volume	0.55 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/80399>

Spec. Version B