

5 µL, Model 75 RN Agilent SYR, Small Removable NDL, 23s ga, 1.71 in, point style AS

Specification Sheet (Part/REF # 87957)



Product Specifications

Volume	5 µL
Syringe Type	Microliter
Syringe Series	700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	23s gauge
Point Style	AS: Conical, Autosampler
Needle Length	1.71 inch (43.4 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.014 inches (0.343 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Needle Inner Diameter	0.010 inches (0.25 mm)
Needle Outer Diameter	0.025 inches (0.64 mm)
Std Needle Dead Volume	0.46 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/87957>

Spec. Version B