

10 µL, Model 701 RN Varian 8100/8200 SYR, Small Removable NDL, 26s ga, 57 mm, point style 5

Specification Sheet (Part/REF # 202880)



Product Specifications

Volume	10 µL
Syringe Type	Microliter
Syringe Series	700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	26s gauge
Point Style	5: Conical Side Port
Needle Length	57 mm
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Needle Inner Diameter	0.005 inches (0.13 mm)
Needle Outer Diameter	0.019 inches (0.474 mm)
Std Needle Dead Volume	0.64 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/202880>

Spec. Version B