

# 10 µL, Model 701 RN Shimadzu SYR, Small Removable NDL, 22s ga, 1.69 in, point style 2

Specification Sheet (Part/REF # 202640)



## Product Specifications

|                        |                         |
|------------------------|-------------------------|
| Volume                 | 10 µL                   |
| Syringe Type           | Microliter              |
| Syringe Series         | 700                     |
| Termination            | Removable Needle (RN)   |
| Needle                 | Small RN                |
| Gauge                  | 22s gauge               |
| Point Style            | 2: Curved, Beveled      |
| Needle Length          | 1.69 in (42.9 mm)       |
| Fluid Path             | SST, borosilicate glass |
| Autoclavable           | No                      |
| Gas Sterilizable       | Yes                     |
| Pressure               | 2000 psig (137.9 bar)   |
| Maximum Temperature    | 50°C / 122°F            |
| Minimum Temperature    | 20°C / 68°F             |
| Needle Inner Diameter  | 0.005 inches (0.13 mm)  |
| Needle Outer Diameter  | 0.028 inches (0.711 mm) |
| Std Needle Dead Volume | 0.76 µL                 |
| Sales Quantity         | 1 syringe per pack      |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/202640>

Spec. Version B