

# 1 mL, Model 1001 LTN CTC SYR (7.6 mm), C-Line, Fixed NDL, 22 ga, point style 3

Specification Sheet (Part/REF # 203081)



## Product Specifications

|                        |                                |
|------------------------|--------------------------------|
| Volume                 | 1 mL                           |
| Syringe Type           | Gastight                       |
| Syringe Series         | 1000                           |
| Termination            | Luer Tip Cemented Needle (LTN) |
| Needle                 | Cemented Needle                |
| Gauge                  | 22 gauge                       |
| Point Style            | 3: Blunt                       |
| Needle Length          | 2.00 in (50.8 mm)              |
| CTC Syringe Type       | C-Line                         |
| Fluid Path             | SST, Borosilicate glass, PTFE  |
| Autoclavable           | No                             |
| Gas Sterilizable       | Yes                            |
| Pressure               | 200 psig (13.8 bar)            |
| Maximum Temperature    | 50°C / 122°F                   |
| Minimum Temperature    | 20°C / 68°F                    |
| Needle Inner Diameter  | 0.0161 inches (0.41 mm)        |
| Needle Outer Diameter  | 0.0283 inches (0.72 mm)        |
| Std Needle Dead Volume | 6.81 µL                        |
| Sales Quantity         | 1 syringe per pack             |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/203081>

Spec. Version B