

10 µL, Model 701 FN CTC SYR (6.6 mm), C-Line, Cemented NDL, 22s ga, point style 3

Specification Sheet (Part/REF # 203073)



Product Specifications

Volume	10 µL
Syringe Type	Microliter
Syringe Series	700
Termination	Fixed Needle (FN)
Needle	Cemented Needle
Gauge	22s gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
CTC Syringe Type	C-Line
Fluid Path	SST, borosilicate glass
Autoclavable	No
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Needle Inner Diameter	0.005 inches (0.13 mm)
Needle Outer Diameter	0.028 inches (0.711 mm)
Std Needle Dead Volume	1.13 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/203073>

Spec. Version B