

100 µL, Model 1710 FN CTC SYR (6.6 mm), C-Line, Fixed NDL, 22s ga, point style 3

Specification Sheet (Part/REF # 203077)



Product Specifications

Volume	100 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Fixed Needle (FN)
Needle	Fixed Needle, no glue or cement
Gauge	22s gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
CTC Syringe Type	C-Line
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Needle Inner Diameter	0.005 inches (0.13 mm)
Needle Outer Diameter	0.028 inches (0.711 mm)
Std Needle Dead Volume	1.13 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/203077>

Spec. Version B