

100 µL, Model 1710 RN CTC SYR, C-Line, NDL Sold Separately

Specification Sheet (Part/REF # 202668)



Product Specifications

Volume	100 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN (CTC)
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
CTC Syringe Type	C-Line
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Std Needle Dead Volume	N/A
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/202668>

Spec. Version B