


**100 µL Gastight CTC PAL
Autosampler Syringe Model 1710 FN
(6.6 mm), X-Type, Fixed Needle, 22
gauge, 2 in, point style 3**

Specification Sheet (Part/REF # 204400)



Product Specifications

Volume	100 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Fixed Needle (FN)
Needle	Fixed Needle, no glue or cement
Gauge	22 gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
CTC Syringe Type	X-Type
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std Needle Dead Volume	6.81 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/204400>

Spec. Version B