

50 µL, Model 1705 N CTC SYR (6.6 mm), S-Line, Cemented Needle, 22 ga, point style 3

Specification Sheet (Part/REF # 67450-01)



Product Specifications

Volume	50 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Cemented Needle (N)
Needle	Cemented Needle
Gauge	22 gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
CTC Syringe Type	S-Line
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.028 inches (0.711 mm)
Std Needle Dead Volume	6.81 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/67450-01>

Spec. Version B