

500 µL, Model 1750 FN CTC SYR (7.9 mm), C-Line, Fixed NDL, 22 ga, point style 3

Specification Sheet (Part/REF # 203349)



Product Specifications

Volume	500 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Fixed Needle (FN)
Needle	Fixed Needle, no glue or cement
Gauge	22 gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
CTC Syringe Type	C-Line
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	700 psig (48.3 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std Needle Dead Volume	6.81 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/203349>

Spec. Version B