

25 µL, Model 1702 CX WISP SYR, 1/4-28 Threads

Specification Sheet (Part/REF # 80020)



Product Specifications

Needle Sold Separately	Yes
Volume	25 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	ChemSeal (C, 1/4"-28)
Needle	No Needle Available
Plunger Style	Waters WISP Plunger
Fluid Path	Borosilicate glass, PTFE, PCTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.029 inches (0.729 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/80020>

Spec. Version B