

10, 25, 50, 100, and 250 µL, Model
RN Waters U6K SYR, 5 Pack,
Removable NDL, 25s ga, 1.97 in,
point style 3

Specification Sheet (Part/REF # 84891)



Syringes for U6K valve come with a 25s gauge needle specifically designed for this valve.

Product Specifications

Volume	10, 25, 50, 100, and 250 µL
Syringe Type	Microliter
Syringe Series	800
Termination	Removable Needle (RN)
Needle	Small and Large RN
Gauge	25s gauge
Point Style	3: Blunt
Needle Length	1.97 in (50.0 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Std Needle Dead Volume	0.91 µL
Sales Quantity	5 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/84891>

Spec. Version B