

1 mL, Model 1001 DAD SYR, Diluter Syringe

Specification Sheet (Part/REF # 201700)



Product Specifications

Needle Sold Separately	Yes
Volume	1 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Diluter AccuDil® (DAD)
Needle	No Needle Available
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/201700>

Spec. Version B