

# 2.5 mL, Model 1002 AD SYR, M8 Threads

Specification Sheet (Part/REF # 201250)



## Product Specifications

|                        |                          |
|------------------------|--------------------------|
| Needle Sold Separately | Yes                      |
| Volume                 | 2.5 mL                   |
| Syringe Type           | Gastight                 |
| Syringe Series         | 1000                     |
| Termination            | AccuDil® (AD)            |
| Needle                 | No Needle Available      |
| Fluid Path             | Borosilicate glass, PTFE |
| Autoclavable           | No                       |
| Gas Sterilizable       | Yes                      |
| Pressure               | 200 psig (13.8 bar)      |
| Maximum Temperature    | 50°C / 122°F             |
| Minimum Temperature    | 20°C / 68°F              |
| Std Needle Dead Volume | N/A                      |
| Sales Quantity         | 1 syringe per pack       |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/201250>

Spec. Version B