

25 mL, Model 1025 AD Syringe, M8 Threads

Specification Sheet (Part/REF # 201400)

The 1000 series syringes are a mid-volume Gastight syringe. This series of syringes is ideal for dispensing volumes from 1 mL up to 100 mL. A variety of different terminations and needle options are available in this series.

Product Specifications

Needle Sold Separately	Yes
Volume	25 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	AccuDil® (AD)
Needle	No Needle Available
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/201400>

Spec. Version B