

250 µL, Model 1725 DAD Syringe, Diluter Syringe

Specification Sheet (Part/REF # 201600)



The 1700 series syringes are the gas tight version of the original 700 series syringes. This series of gas tight syringes is ideal for dispensing volumes from 1 µL up to 500 µL. A variety of different terminations and needle options are available in this series.

Product Specifications

Needle Sold Separately	Yes
Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Diluter AccuDil® (DAD)
Needle	No Needle Available
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Std Needle Dead Volume	N/A
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/201600>

Spec. Version B