

500 µL, Model 1750 AD Syringe, M8 Threads

Specification Sheet (Part/REF # 201150)



The 1700 series syringes are the gas tight version of the original 700 series syringes. This series of gas tight syringes is ideal for dispensing volumes from 1 µL up to 500 µL. A variety of different terminations and needle options are available in this series.

Product Specifications

Needle Sold Separately	Yes
Volume	500 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	AccuDil® (AD)
Needle	No Needle Available
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	700 psig (48.3 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Std Needle Dead Volume	N/A
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/201150>

Spec. Version B