

# 500 mL, Model S0500 Tracheal Syringe, No Needle Available

Specification Sheet (Part/REF # 86301)



Super Syringes are extremely high volume Gastight syringes. They work best with volumes between 500 mL to 2 L.

## Product Specifications

<b>Needle Sold Separately</b>	Yes
<b>Volume</b>	500 mL
<b>Syringe Type</b>	Gastight
<b>Syringe Series</b>	Super Syringe
<b>Termination</b>	Tracheal
<b>Needle</b>	No Needle Available
<b>Gauge</b>	No Needle Included
<b>Point Style</b>	No Needle Included
<b>Needle Length</b>	No Needle Included
<b>Fluid Path</b>	Buna-N, Nickel plated brass, PTFE, Silicone, Acrylic
<b>Autoclavable</b>	No
<b>Gas Sterilizable</b>	Yes
<b>Pressure</b>	30 psig (2.1 bar)
<b>Maximum Temperature</b>	80°C / 176°F
<b>Minimum Temperature</b>	0°C / 32°F
<b>Barrel Inner Diameter</b>	2.44 inches (62.0 mm)
<b>Barrel Outer Diameter</b>	2.95 inches (74.9 mm)
<b>Std Needle Dead Volume</b>	N/A
<b>Sales Quantity</b>	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/86301>

**Spec. Version B**