

10, 25, 50, 100, and 250 µL, Model RN SYR, 5 Pack, Removable NDL, 26s and 22s ga, 2 in, point style 2

Specification Sheet (Part/REF # 84890)



The 800 series syringes are the reinforced plunger version of the original 700 series syringes. Attached to the flange of the syringe is a blue Syringe Holder. The Syringe Holder provides a spot to hold the syringe that will not heat up the barrel causing dispense inaccuracies.

Product Specifications

Volume	10, 25, 50, 100, and 250 µL
Syringe Type	Microliter
Syringe Series	800
Termination	Removable Needle (RN)
Needle	Small and Large RN
Gauge	26s and 22s gauge
Point Style	2: Curved, Beveled
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Std Needle Dead Volume	N/A
Sales Quantity	5 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/84890>

Spec. Version B