

250 µL, Model 1725 DX, Diluter Syringe

Specification Sheet (Part/REF # 81126)



Product Specifications

Needle Sold Separately	Yes
Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Diluter PTFE Luer Lock (D)
Needle	No Needle Available
Plunger Style	X-Style
Fluid Path	Borosilicate glass, PTFE, PCTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/81126>

Spec. Version B