

500 µL, Model 1750 DX, Diluter Syringe

Specification Sheet (Part/REF # 81226)



Product Specifications

Needle Sold Separately	Yes
Volume	500 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Diluter PTFE Luer Lock (D)
Needle	No Needle Available
Plunger Style	X-Style
Fluid Path	Borosilicate glass, PTFE, PCTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	700 psig (48.3 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.128 inches (3.26 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/81226>

Spec. Version B