

10 mL, Model 1010 Bubble Free Prime (BFP), Instrument Syringe

Specification Sheet (Part/REF # 59000-50)



The Bubble Free Prime syringes are compatible with all Microlab 600 applications. The SaltLine syringes were specially developed to withstand harsh working conditions in demanding applications over a long product lifetime.

Product Specifications

| | |
|------------------------------|------------------------------------|
| Volume | 10 mL |
| Syringe Type | Gastight, Bubble Free Prime |
| Syringe Series | 1000 |
| Termination | Bubble Free Prime™ (BFP, 5/16"-28) |
| Needle | No Needle Available |
| Plunger Style | X-Style |
| Fluid Path | Borosilicate glass, PTFE, PFA |
| Autoclavable | No |
| Gas Sterilizable | Yes |
| Pressure | 200 psig (13.8 bar) |
| Maximum Temperature | 80°C / 176°F |
| Minimum Temperature | 10°C / 50°F |
| Barrel Inner Diameter | 0.574 inches (14.6 mm) |
| Barrel Outer Diameter | 0.695 inches (17.7 mm) |
| Sales Quantity | 1 syringe per pack |
| Plunger Tip Material | PTFE |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/59000-50>

Spec. Version B