

10 µL, Model 801 RN Waters U6K SYR, Small Removable NDL, 25s ga, 1.97 in, point style 3

Specification Sheet (Part/REF # 84815)



Syringes for U6K valve come with a 25s gauge needle specifically designed for this valve.

Product Specifications

Volume	10 µL
Syringe Type	Microliter
Syringe Series	800
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	25s gauge
Point Style	3: Blunt
Needle Length	1.97 in (50.0 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.019 inches (0.485 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Needle Inner Diameter	0.007 Inches (0.18 mm)
Needle Outer Diameter	0.02 inches (0.508 mm)
Std Needle Dead Volume	0.91 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/84815>

Spec. Version B