

100 µL, Model 1710 TLL SYR, Long Life Instrument Syringe

Specification Sheet (Part/REF # 9882-20)



Product Specifications

Volume	100 µL
Syringe Type	Longlife Gastight
Syringe Series	1700
Termination	PTFE Luer Lock (TLL)
Plunger Style	Side Post Mounting (XL)
Fluid Path	Borosilicate glass, PTFE, UHMWPE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/9882-20>

Spec. Version B