

# 250 µL, Model 1725 TLL SYR, Long Life Instrument Syringe

Specification Sheet (Part/REF # 8300-25)



## Product Specifications

|                       |                                  |
|-----------------------|----------------------------------|
| Volume                | 250 µL                           |
| Syringe Type          | Longlife Gastight                |
| Syringe Series        | 1700                             |
| Termination           | PTFE Luer Lock (TLL)             |
| Plunger Style         | X-Style                          |
| Fluid Path            | Borosilicate glass, PTFE, UHMWPE |
| Autoclavable          | No                               |
| Gas Sterilizable      | Yes                              |
| Pressure              | 500 psig (34.5 bar)              |
| Maximum Temperature   | 80°C / 176°F                     |
| Minimum Temperature   | 10°C / 50°F                      |
| Barrel Inner Diameter | 0.091 inches (2.30 mm)           |
| Barrel Outer Diameter | 0.305 inches (7.75 mm)           |
| Sales Quantity        | 1 syringe per pack               |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/8300-25>

Spec. Version B