

# 250 µL, Model 1725 TLL SYR,

## Instrument Syringe

Specification Sheet (Part/REF # 9880-25)



Designing a system that needs a syringe? Utilize our syringe experts to get the best connections and materials of construction for your application.

### Product Specifications

Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	PTFE Luer Lock (TLL)
Needle	Instrument Syringe
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Plunger Style	Side Post Mounting (XL)
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Std Needle Dead Volume	N/A
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/9880-25>

Spec. Version B