

10 mL, Model 1010 LTSN Syringe, Cemented Needle, Custom gauge, length, point style

Specification Sheet (Part/REF # 81614)



The 1000 series syringes are a mid-volume Gastight syringe. This series of syringes is ideal for dispensing volumes from 1 mL up to 100 mL. A variety of different terminations and needle options are available in this series.

Product Specifications

| | |
|------------------------|---|
| Volume | 10 mL |
| Syringe Type | Gastight |
| Syringe Series | 1000 |
| Termination | Luer Tip Cemented Needle Special (LTSN) |
| Needle | Cemented Needle |
| Gauge | Custom Gauge |
| Point Style | 2, 3, 4, 5, or AS |
| Needle Length | Custom Needle Length |
| Fluid Path | SST, Borosilicate glass, PTFE |
| Autoclavable | No |
| Gas Sterilizable | Yes |
| Pressure | 200 psig (13.8 bar) |
| Maximum Temperature | 50°C / 122°F |
| Minimum Temperature | 10°C / 50°F |
| Barrel Inner Diameter | 0.574 inches (14.6 mm) |
| Barrel Outer Diameter | 0.695 inches (17.7 mm) |
| Std Needle Dead Volume | N/A |
| Sales Quantity | 1 syringe per pack |

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/81614>

Spec. Version B