

Custom Length PTFE Tubing, 11 Gauge (0.120" O.D. / 0.096" I.D.)

Specification Sheet (Part/REF # 20911)



Product Specifications

| | |
|------------------------------|---------------------|
| Tubing Material | PTFE |
| Gauge | 11 gauge |
| Tubing Inner Diameter (I.D.) | 0.096 in / 2.438 mm |
| Tubing Outer Diameter (O.D.) | 0.120 in / 3.048 mm |
| Tubing Length | Custom Length |
| Fitting Material | No Hub |
| Sales Quantity | 1 each |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/fittings/20911>

Spec. Version B