

Custom Length PTFE Tubing, 13 Gauge (0.101" O.D. / 0.077" I.D.)

Specification Sheet (Part/REF # 20913)



Product Specifications

Tubing Material	PTFE
Gauge	13 gauge
Tubing Inner Diameter (I.D.)	0.077 in / 1.956 mm
Tubing Outer Diameter (O.D.)	0.101 in / 2.565 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/fittings/20913>

Spec. Version B