

Custom Length PTFE Tubing, 14 Gauge (0.092" O.D. / 0.068" I.D.)

Specification Sheet (Part/REF # 20914)



Product Specifications

Tubing Material	PTFE
Gauge	14 gauge
Tubing Inner Diameter (I.D.)	0.068 in / 1.727 mm
Tubing Outer Diameter (O.D.)	0.092 in / 2.337 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/fittings/20914>

Spec. Version B