

# Custom Length PTFE Tubing, 15 Gauge (0.085" O.D. / 0.061" I.D.)

Specification Sheet (Part/REF # 20915)



## Product Specifications

Tubing Material	PTFE
Gauge	15 gauge
Tubing Inner Diameter (I.D.)	0.061 in / 1.549 mm
Tubing Outer Diameter (O.D.)	0.085 in / 2.159 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/fittings/20915>

Spec. Version B