

Custom Length PTFE Tubing, 26 Gauge (0.037" O.D. / 0.019" I.D.)

Specification Sheet (Part/REF # 20926)



Product Specifications

Tubing Material	PTFE
Gauge	26 gauge
Tubing Inner Diameter (I.D.)	0.019 in / 0.483 mm
Tubing Outer Diameter (O.D.)	0.037 in / 0.940 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/fittings/20926>

Spec. Version B