

Custom Length PTFE Tubing, 8 Gauge (0.164" O.D. / 0.134" I.D.)

Specification Sheet (Part/REF # 20908)



Product Specifications

Tubing Material	PTFE
Gauge	8 gauge
Tubing Inner Diameter (I.D.)	0.134 in / 3.404 mm
Tubing Outer Diameter (O.D.)	0.164 in / 4.166 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/fittings/20908>

Spec. Version B