

**Female Luer Lock, One Fitting with  
Custom Length PTFE Tubing, 14  
Gauge (0.092" O.D. / 0.068" I.D.)**

Specification Sheet (Part/REF # 90664)



**Product Specifications**

|                              |                     |
|------------------------------|---------------------|
| Fitting Material             | CTFE                |
| Tubing Material              | PTFE                |
| Gauge                        | 14 gauge            |
| Tubing Inner Diameter (I.D.) | 0.068 in / 1.727 mm |
| Tubing Outer Diameter (O.D.) | 0.092 in / 2.337 mm |
| Tubing Length                | Custom Length       |
| Number of Fittings           | 1                   |
| Fitting Type                 | Female Luer         |
| Sales Quantity               | 1 each              |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/fittings/90664>

**Spec. Version B**