

PCR FramePlate® 384-well

Specification Sheet (Part/REF # 814305)



Hamilton 384-well PCR FramePlate®

Product Specifications

Maximum Well Capacity 50 µL	
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Barcode Label	Yes
Frame Material	Polycarbonate - PC
Tube Material	Polypropylene - PP
Recycling Code	7
SLAS ANSI	Yes
Packaging	10 plates per bag
Barcode	No
Compatible Carrier	SLAS/ANSI carrier
Quantity	50 plates per case
Disclaimer	PCR FramePlate® is covered by one or more of the following US patents or their foreign counterparts, owned by Eppendorf AG: US Patent Nos. 7,347, 977 and 6,340,589. PCR FramePlate® is a registered trademark owned by 4titude® Ltd.

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/814305>

Spec. Version B