

1.2 mL Deep Well Plate - Only Available in the U.S.

Specification Sheet (Part/REF # 6472-01)



1.2 mL Deep Well Plate; Polystyrene

Product Specifications

Maximum Well Capacity	1.2 mL
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, ATP, and PCR inhibitors
Material	Polystyrene - PS
Recycling Code	6
SLAS ANSI	Yes
Packaging	4 plates per bag
Barcode	No
Compatible Carrier	SLAS/ANSI carrier
Quantity	32 plates per case

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/1669891>

Spec. Version B