


Hamilton PCR ComfortLid

Specification Sheet (Part/REF # 814300)



Hamilton PCR two-component ComfortLid

Product Specifications

Purity Level	Sterile, Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, PCR inhibitors
Material	PC - Polycarbonate, TPE - Thermoplastic Elastomers
Packaging	10 lids per blister
Barcode	No
Quantity Per Case	50 lids per case
Compatible Carrier	814337 PCR ComfortLid Adapter for blister on deck integration

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/814300>

Spec. Version B