

1.5 mL FlipTube Standard

Specification Sheet (Part/REF # 235454)



FlipTube Standard Non-sterile

Product Specifications

Platform Compatibility	Microlab VANTAGE, Microlab STAR
Volume	1.5mL
Sterile	Non-Sterile
Color	Clear
Material	Polypropylene - PP
Recycling Code	5
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Temperature Range	From -80 °C (max. 1 mL) to +55 °C (max 1.5 mL)
Tubes / Case	2000 pcs

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/235454>

Spec. Version B