

1.5 mL FlipTube Biological Purityplus Sterile

Specification Sheet (Part/REF # 235693)



1.5 mL FlipTube Biological Purityplus, Sterile

Product Specifications

Platform Compatibility	Microlab VANTAGE, Microlab STAR
Volume	1.5 mL
Sterile	Sterile
Color	Clear
Material	Polypropylene - PP
Recycling Code	5
Purity Level	Sterile Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, PCR inhibitors, and ATP
Temperature Range	From -80 °C (max. 1 mL) to +55° C (max 1.5 mL)
Tubes/ Case	2000 pcs

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/other-robotics/235693>

Spec. Version B