

60 mL Black Reagent Reservoir, Self-Standing with Lid, Bulk

Specification Sheet (Part/REF # 56694-03)



60 mL Black Reagent Reservoir, Self-Standing with Lid, No Barcode

Product Specifications

Material	Polypropylene - PP
Production Compliance	Product produced in a Class 8 Clean Room
Purity Level	Free of DNase, RNase, Human DNA, and PCR inhibitors
Recycling Code	5
Quantity Per Case	Case of 20 reagent reservoirs
Packaging	20 reservoirs per bag
SLAS ANSI	No
Barcode	No
Compatible Carrier	53646-01 60 mL Reagent Reservoir Carrier - 5 reservoir positions

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/56694-03>

Spec. Version B