

50 mL Reagent Reservoir

Specification Sheet (Part/REF # 187297)



50 mL Reagent Reservoir Non-Self-Standing without Lid, No Barcode

Product Specifications

Material	Polypropylene - PP
Purity Level	Standard Production
Recycling Code	5
Quantity Per Case	Case of 12 reservoirs
Barcode	No
Compatible Carrier	188047 MultiFlex 50 mL Reagent Trough Module - 6 reservoir positions, 187299 50 mL Reagent Reservoir Carrier - 5 reservoir positions

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/187297>

Spec. Version B