

120 mL Reagent Reservoir Self- Standing with Barcode Label

Specification Sheet (Part/REF # 194052)



120 mL Reagent Reservoir Self-Standing with Lid, Barcode Label included

Product Specifications

Material	Polypropylene - PP
Production Compliance	Product produced in a Class 8 Clean Room
Purity Level	Free of DNase, RNase
Recycling Code	5
Quantity per Case	Case of 12 reagent reservoirs
Packaging	Individually wrapped
Barcode	Yes
Compatible Carrier	194058 120 mL Reagent Reservoir Carrier - 3 reservoir positions

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/194502>

Spec. Version B