

120 mL Reagent Reservoir with V-Bottom

Specification Sheet (Part/REF # 182703)



120 mL Reagent Reservoir with V-Bottom, Non-Self-Standing, Barcode Label included without Lid

Product Specifications

Material	Polypropylene - PP
Purity	Standard production
Recycling Code	5
Quantity per Case	Case of 12 reagent reservoirs
Packaging	Individually wrapped
Barcode	No
SLAS ANSI	No
Compatible Carrier	185290 120 mL Reagent Reservoir Carrier - 1T holds 3x 120 mL reservoirs carrier

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/182703>

Spec. Version B